

*Document Detail*

9/25/2019

1

**Type:** INSERT  
**Document No.:** 3IX762[08]  
**Title:** HBT (HETEROPHILIC BLOCKING TUBE)  
**Owner:** HRODRIGUEZ Homar Rodriguez Gonzalez  
**Status:** CURRENT  
**Effective Date:** 16-Nov-2010 **Expiration Date:** 16-Nov-2023

Revision Notes

<u>Document Build No.</u>	<u>Access Activity</u>	<u>Accessed By</u>	<u>Accessed Date</u>
1	Check In	JSOMERS	15-Nov-2010

Note:

Review

**Build No.:** 1 **Closed Date:** 11/16/2010 4:37:01PM

**Review:** Change Order Document Review - Auto Close

**Review Purpose:** Take document from Inworks status to Current status with Auto Close option

**Review Note:** SYSTEM AUTO CLOSE REVIEW

<u>Level</u>	<u>Owner Role</u>	<u>Actor</u>	<u>Sign-off Date</u>	<u>Sign-off By</u>
1	CO Approver - ICD CO Approver - ICD	HTA Hien Ta	15-Nov-2010 10:40 pm	HTA

**Note To Approver:**

**Note From Approver:**Approved

2	Document Owner Document Owner	TTRAN Tuan Tran	16-Nov-2010 12:01 am	TTRAN
---	-------------------------------	-----------------	----------------------	-------

**Note To Approver:**

**Note From Approver:**Approved

3	CO Approver - QA CO Approver - QA	RBEARRY Roy Bearry	16-Nov-2010 12:04 am	RBEARRY
---	-----------------------------------	--------------------	----------------------	---------

**Note To Approver:**

**Note From Approver:**Reviewed and Approved

4	DOC CONTROL Document Control	DOCCONTROL Document Control	16-Nov-2010 4:37 pm	JSOMERS
---	------------------------------	-----------------------------	---------------------	---------

**Note To Approver:**

**Note From Approver:**Reviewed and Approved

# HBT (Heterophilic Blocking Tube)

# Scantibodies

Elimination of Heterophilic Interference for Immunochemistry Assays

(Part Number: 3IX762)


For Investigational Use Only

Store at 2° - 8° C


## THE HETEROPHILIC INTERFERENCE PROBLEM

A heterophilic sample is a serum or plasma sample which contains antibodies which are able to bind to animal antibodies used in immunochemistry assays. The most commonly reported assay interference effect of heterophilic antibodies is a false positive assay result. False negative assay results have also been reported in the literature.

The following diagram illustrates a normal sandwich immunoassay where the concentration of the analyte is responsible for the positive assay result.



The following diagram illustrates a sandwich assay where the heterophilic antibody is responsible for the false positive assay result.



It has been found that as much as 22% of certain sandwich immunoassay results are false positive results caused by heterophilic antibody interference. With such a large potential for immunoassay false positive values it is important to block heterophilic interference.

## INTENDED USE

The HBT allows for the rapid and simple elimination of false positive heterophilic interference in plasma or serum for immunoassays (i.e., FSH, LH, Prolactin, TSH, Ferritin, CEA, AFP, hCG, HBsAg, CK-MB, CA 125<sup>1</sup>, CA 19-9<sup>1</sup>, NSE, etc.). HBT represents a sample pretreatment/second assay intended to confirm or disqualify the original FDA licensed non pretreatment assay result. The assay result from the pretreatment is **NEVER TO BE USED AS A REPORTABLE RESULT**. The pretreatment is only a confirmation aid designed to assist the lab to know whether to report the original non pretreatment assay result. In other words, if the pretreatment assay result is the same as the original assay result the original result is reported. However, if the pretreatment sample assay result is lower than the original result, the original result is not reported and the sample is submitted for further study for potential false positive assay interference.

## REAGENT

The HBT contains a unique blocking reagent composed of specific binders which inactivate heterophilic antibodies. Once the specific binders have bound to the heterophilic antibodies, the antibodies are no longer able to cause immunoassay interference. Each tube contains enough reagent to inactivate the heterophilic antibodies in 500  $\mu$ L of sample. The reagent is in the form of a lyophilized pellet at the bottom of the tube.

## STORAGE AND STABILITY

Upon receipt, store the tubes at 2° - 8° Celsius.

**NOTE:** Studies have shown that there is no significant impact on HBT performance when temporarily storing the tubes for up to 6 weeks at room temperature or below - 20° C.

## MATERIALS REQUIRED BUT NOT PROVIDED

Pipetting device with tips for 500  $\mu$ L

<sup>1</sup>CA 19-9 and CA 125 are trademarks of Centocor, Inc.

## PROCEDURE

1. Use one tube for each sample.
2. Holding the HBT upright, gently tap the bottom of the tube on a hard surface. This action brings all of the reagent to the bottom of the tube.
3. Remove the cap from the tube.
4. Pipette 500  $\mu$ L of the patient sample into the bottom of the tube.
5. Avoid sample carryover contamination by using a new pipette tip for each sample.
6. Cap the tube and invert 5 times to mix the sample with the reagent.
7. Incubate for 1 hour at room temperature (18° - 28° C).
8. Assay the now treated sample to obtain a result that is free from heterophilic antibody interference.

Note: When testing samples of plasma, some samples may develop a gel like consistency when first mixed in the tube. In most cases, gentle agitation of the tube will liquefy the sample so that it may be tested.

## PRECAUTIONS

1. For investigational use only.
2. Once a patient sample has been added to the HBT, do not interchange caps among different tubes.
3. A new dispensing tube or pipette tip should be used for each patient sample.

## INTERPRETATION OF RESULTS

If the HBT is used for the initial assay, consider the immunochemistry result to be unaffected by heterophilic antibody interference.

If the HBT is used for a secondary confirmation assay, compare the results from the first assay (initial sample not treated with HBT) and the confirmation assay (second sample treated with HBT). If the assay result from the HBT treated sample is different from the assay result from the untreated sample, the difference is due to heterophilic interference.

## LIMITATIONS OF PROCEDURE

1. For diagnostic purposes, the results obtained by this sample treatment should be used as an adjunct to other data (e.g., symptoms, results of other testing, clinical impression, etc.) available to the physician.
2. There may be some samples with extremely strong heterophilic interference. In such cases the HBT may not be able to block all of the assay interference.
3. The HBT is designed for application in analyte determinations.

**SCANTIBODIES LABORATORY, INC.**  
9336 Abraham Way, Santee, CA 92071  
Phone 1-619-258-9300, Fax 1-619-258-9366