

PRODUCT SPECIFICATION SHEET

DESCRIPTION: **SYSCON TUMOR MARKERS 1**

PRODUCT NUMBER: **3CA064**

| ANALYTE | UNIT | ANALYTE RANGE | ASSAY METHOD |
|----------------------|-------------|----------------------|---------------------|
| CEA | ng/mL | 4 - 6 | Elecsys |
| AFP | IU/mL | 6 - 10 | Elecsys |
| CA 19-9 | IU/mL | 15 - 30 | Elecsys |
| CA 125 | IU/mL | 25 - 35 | Elecsys |
| CA 15-3 | IU/mL | 15 - 30 | Elecsys |
| CA 72-4 | IU/mL | 4 - 6 | Elecsys |
| CYFRA | ng/mL | 2 - 4 | Elecsys |
| Ferritin | ng/mL | 20 - 30 | Elecsys |
| β-HCG | mIU/mL | 5 - 20 | Elecsys |
| NSE | μg/L | 10 - 15 | Elecsys |
| Thyroglobulin | ng/mL | 30 - 50 | Elecsys |
| PSA | ng/mL | 3 - 5 | Elecsys |
| Free PSA | ng/mL | 0.5 - 2 | Elecsys |

MICROBIOLOGY: ≤ 100 cfu/mL (Colony Forming Units per mL)

STORAGE: 2-8°C

SHELF LIFE: 24 months from date of manufacture

PRODUCT SPECIFICATION SHEET

DESCRIPTION: **SYSCON TUMOR MARKERS 2**

PRODUCT NUMBER: **3CB064**

| ANALYTE | UNIT | ANALYTE RANGE | ASSAY METHOD |
|----------------------|-------------|----------------------|---------------------|
| CEA | ng/mL | 40 - 60 | Elecsys |
| AFP | IU/mL | 50 - 150 | Elecsys |
| CA 19-9 | IU/mL | 80 -120 | Elecsys |
| CA 125 | IU/mL | 80 - 120 | Elecsys |
| CA 15-3 | IU/mL | 80 - 120 | Elecsys |
| CA 72-4 | IU/mL | 30 - 70 | Elecsys |
| CYFRA | ng/mL | 20 - 40 | Elecsys |
| Ferritin | ng/mL | 160 - 240 | Elecsys |
| β-HCG | mIU/mL | 500 -1500 | Elecsys |
| NSE | μg/L | 60 - 140 | Elecsys |
| Thyroglobulin | ng/mL | 180 - 220 | Elecsys |
| PSA | ng/mL | 30 - 50 | Elecsys |
| Free PSA | ng/mL | 5 - 15 | Elecsys |

MICROBIOLOGY: ≤ 100 cfu/mL (Colony Forming Units per mL)

STORAGE: 2-8°C

SHELF LIFE: 24 months from date of manufacture