

PRODUCT SPECIFICATION SHEET

DESCRIPTION: **REDICON SP LOW LEVEL**

PRODUCT NUMBER: **3CA416**

ANALYTE	UNIT	ANALYTE RANGE	ASSAY METHOD
CEA	ng/mL	2.5 - 3.0	AXSYM
AFP	ng/mL	27 - 33	AXSYM
IgG	mg/dL	720- 800	ARRAY
IgA	mg/dL	72 - 88	ARRAY
IgM	mg/dL	31.5 - 38.5	ARRAY
C3	mg/dL	72 - 88	ARRAY
C4	mg/dL	15.3 - 18.7	ARRAY
CRP	mg/dL	0.63 - 0.77	ARRAY
RF	IU/mL	22.5 - 27.5	ARRAY
ASO	IU/mL	135 - 165	ARRAY

MICROBIOLOGY: ≤ 100 cfu/mL (Colony Forming Units per mL)

STORAGE: 2-8°C

SHELF LIFE: 24 months from date of manufacture

PRODUCT SPECIFICATION SHEET

DESCRIPTION: **REDICON SP HIGH LEVEL**

PRODUCT NUMBER: **3CB416**

ANALYTE	UNIT	ANALYTE RANGE	ASSAY METHOD
CEA	ng/mL	16.2 - 19.8	AXSYM
AFP	ng/mL	162 - 198	AXSYM
IgG	mg/dL	1530 - 1870	ARRAY
IgA	mg/dL	225 - 275	ARRAY
IgM	mg/dL	162 - 198	ARRAY
C3	mg/dL	153 - 187	ARRAY
C4	mg/dL	31.5 - 38.5	ARRAY
CRP	mg/dL	4.5 - 5.5	ARRAY
RF	IU/mL	63 - 77	ARRAY
ASO	IU/mL	243 - 297	ARRAY

MICROBIOLOGY: ≤ 100 cfu/mL (Colony Forming Units per mL)

STORAGE: 2-8°C

SHELF LIFE: 24 months from date of manufacture